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MONTES L. D.; ROCKLEY, M. G. **Teacher Perceptions in the Selection of Experiments**. Journal of Chemical Education, v. 79, n. 2, 244-248. feb. 2002.

RANDALL, J. **Advanced Chemistry with Vernier 4th Edition**. Vernier Software & Technology. *Beaverton, OR*. 2019.

# Electroplating

# 21

In this experiment, you will conduct, observe, and measure the process of electroplating. This process is used to deposit a layer of metal, such as chromium, copper, or gold, onto another metal. As a commercial process, electroplated coatings are used to improve appearance, resist corrosion, or improve hardness of metallic surfaces. This experiment describes one method of producing a copper coating on a brass key or other suitable metallic object.

You will prepare an electrochemical cell by using a copper strip as the cathode (positive terminal) and a brass key as the anode (negative terminal). The electrodes are immersed in a solution containing acidified copper (II) sulfate. As you apply a potential to the electrodes, you will be effectively transferring Cu atoms from the anode to the surface of the brass key.

In this experiment, you will use one application of Faraday's law, stated in equation form below.

$$\text{Mass deposited at an electrode} = \frac{I \times t \times (MM)}{96,500 \text{ C/mol } e^- \times n}$$

$I$  is the current in amperes(A);  $t$  is the time that the current is applied, in seconds;  $MM$  is the molar mass of the element that is deposited;  $n$  is the number of moles of electrons/mol; and  $96,500 \text{ C/mol } e^-$  is  $F$ , the Faraday constant.


Figure 1

## OBJECTIVES

In this experiment, you will

- Prepare and operate an electrochemical cell to plate copper onto a brass surface.
- Measure the amount of copper that was deposited in the electroplating process.
- Calculate the amount of energy used to complete the electroplating process.

## MATERIALS

Vernier computer interface	electrolyte solution (CuSO <sub>4</sub> in H <sub>2</sub> SO <sub>4</sub> )
computer	vinegar
Vernier Constant Current System	1 cm × 10 cm copper strip
or 1.5 volt DC power supply and	brass key
Vernier Current Probe	solid sodium chloride, NaCl
steel wool	balance, 0.001 g precision
two 250 mL beakers	bare copper wire, 20-22 gauge
distilled water	connecting wires with alligator clips
	(for use with the Current Probe only)

## PROCEDURE

1. Obtain and wear goggles.
2. Use steel wool to clean a brass key and a strip of copper, which will be the electrodes of the electrochemical cell.
3. Mix 3 g of NaCl with 15 mL of vinegar in a 250 mL beaker. Wash the key and the copper strip in this salt-vinegar solution. Rinse the key and copper strip with distilled water and dry each metal piece.
4. Use an analytical balance to determine the mass of the key and the mass of the copper strip. Record these two masses in your data table.
5. Fill a 250 mL beaker about 3/4 full with the electrolyte solution. **CAUTION:** *The electrolyte solution in this experiment is prepared in H<sub>2</sub>SO<sub>4</sub> and should be handled with care.*
6. Attach a 7 cm length of bare copper wire to the brass key to act as a handle, so that it is easier to completely immerse the key in the electrolyte solution. Do not immerse until Step 10.
7. Gently turn the dial of the Constant Current System counterclockwise to confirm that it is in the minimum current position. Clip the wire handle on the key (anode) and the copper strip (cathode) to the proper leads (see Figure 1). **Important:** You will not place the electrodes in the beaker until Step 10.